

# **BloodSTOP**<sup>®</sup>

**CONTROLS BLEEDING FAST**


## **BloodSTOP<sup>®</sup> Instructions for Use**

BloodSTOP<sup>®</sup> Hemostatic Wound Dressing is made of 100% natural etherified sodium carboxymethyl cellulose, which is water soluble and biocompatible. When BloodSTOP material comes into contact with blood and exudates, it turns into a hemostatic gel that adheres to the wounded area to promote healing. The clear gel prevents contamination and enables easy wound monitoring.


### **Indications For Use:**

For topical wound bleeding control from lacerations, cuts, abrasions, and venipuncture, and for fast wound healing when a moist wound environment is needed. For the general population and (EU countries only) for those on anticoagulant medications.


### **Directions:**



1. Clean wound. Apply 1 or more layers to wound with dry or gloved hand. Apply continuous pressure until bleeding stops completely, usually 2-3 minutes. To avoid sticking to hand, roll fingers off of sticky gel.



2. Cover with additional dressing to maintain moist environment and pressure. If bleeding continues, seek medical attention.



3. Change within 24 hours. To remove, saturate with sterile water or sterile saline and gently wipe away the gel, being careful not to disturb the clot. Repeat above steps until wound heals.

### **Nose Bleeding:**



1. Insert gently by hand, or with cotton swab, as close to the source of bleeding as possible.

2. Compress both nostrils from above until bleeding stops. Repeat if necessary. Seek medical attention if bleeding does not stop.



3. To remove, gently tug or squeeze gel out of the nostril. Gently rinse with sterile water or sterile saline.

**WARNING:** In case of swelling or redness, discontinue use. For topical use only. Cover BloodSTOP with a wrap or adhesive bandage to maintain pressure if the patient is to be moved. Do not reuse, as device will no longer be sterile after initial use, and could result in infection.

**STORAGE:** For maximum shelf life, storage temperature of 15-30 °C (59-86 °F) is recommended.

**INGREDIENTS:** 100% natural etherified sodium carboxymethyl cellulose.



WellKang Ltd.  
Enterprise Hub, NW Business Complex,  
1 Beraghmore Rd. Derry, BT48 8SE,  
Northern Ireland, UK



*Not made with natural rubber latex*

### Symbols Used on Labeling



Single use



Use until year and month



Reference number



See instructions for use



Batch number



Do not use if package is damaged



Company address




Do not resterilize



Temperature Limits



Sterile unless package is opened or damaged.  
Method of sterilization: irradiation



CE Mark and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC as amended by 2007/47/EEC.



### LifeSciencePLUS, Inc.

2520A Wyandotte Street, Mountain View, CA 94043, U.S.A.

Mailing: P.O. Box 60783, Palo Alto, CA 94306, U.S.A

Tel: +1 (650) 565-8172

[www.lifescienceplus.com](http://www.lifescienceplus.com); [www.bloodstopix.com](http://www.bloodstopix.com)

© 2021 LifeScience PLUS, Inc. All Rights Reserved.

BloodSTOP® is a registered trademark. BloodSTOP® and

BloodSTOP® iX are covered by U.S. Patent No. 7,262,181

Manufactured in the U.S.A.