


**BloodSTOP® iX**  
CONTROLS BLEEDING FAST

**ADVANCED HEMOSTAT  
AND WOUNDHEAL®**

*Power to control bleeding  
Power to heal*



## **BloodSTOP® iX Absorbable Hemostat and WoundHEAL®**

- ▶ Class III CE Mark approval for implant surgical use
- ▶ Hemostasis for patients with anti-coagulant medication therapy
- ▶ 100% natural, plant-based absorbable matrix without any residue



- Controls bleeding fast
- At the wound interface, upon contact with blood, exudate, or other liquid, BloodSTOP® iX matrix quickly turns into a gel and seals the wound
- The gel forms a barrier to seal and protect the wound, and quickly activates the intrinsic coagulation pathway, thus accelerating clot formation
- Minimizes blood loss and provides an ideal environment to promote fast wound healing and tissue growth
- Very cost-effective

**Burn & Trauma  
Surgery**

**Neurosurgery**

**Cardiovascular Surgery**

**General Surgery**

**Thoracic Surgery**

**Gynecologic Surgery**

**Urologic Surgery**

**LifeScience PLUS**  
REVOLUTIONARY WOUND CARE



**BloodSTOP® iX**  
CONTROLS BLEEDING FAST



## ADVANCED HEMOSTAT AND WOUNDHEAL®

When applied on a surgical wound, BloodSTOP® iX quickly absorbs blood and exudates, transforms into a gel to seal the wound with a protective transparent layer, acts fast to control bleeding, and creates an environment for wound healing. The clear gel facilitates monitoring of hemostasis.



BloodSTOP iX Used in Thoracic Surgery



BloodSTOP iX Used in Intracranial Surgery



BloodSTOP iX Used in Laparoscopic Cyst Reduction Surgery

BloodSTOP® iX Absorbable Hemostat and WoundHEAL® is used for surgical procedures to quickly control capillary, venous, and small arterial hemorrhage. It can be used in neurosurgery, cardiovascular surgery, general surgery, thoracic surgery, gynecologic surgery, urologic surgery, and burn & trauma surgery.

### Features & Benefit

- Controls bleeding fast, turns into gel to seal the wound
- Accelerates clot formation
- Translucent gel – easy to monitor
- Creates moist environment to promote healing
- No animal, foreign human, or added-chemical components
- Easy to cut and apply, conforms to any wound
- pH neutral
- 100% absorbable without any residue
- 100% natural, biocompatible, non-irritating
- Very cost effective

Product Description	Size	Quantity Per Box	Product Code
BloodSTOP® iX Absorbable Hemostat and WoundHEAL®	0.5"x2" (1.3x5cm)	24/box	BS-iX-27
	2"x2" (5x5cm)	12/box	BS-iX-14
	2"x4" (5x10cm)	12/box	BS-iX-15
	4"x8" (10x20cm)	12/box	BS-iX-17
	2"x14" (5x35cm)	12/box	BS-iX-20

#### LifeSciencePLUS, Inc.

2520A Wyandotte Street, Mountain View, CA 94043, USA

Mailing: P.O. Box 60783, Palo Alto, CA 94306

1.877.587.5433 (U.S.) 1.650.565.8172 (International)

www.lifescienceplus.com

www.bloodstopix.com

© 2019 LifeScience PLUS, Inc. All Rights Reserved.

BloodSTOP® is a registered trademark.

BloodSTOP® and BloodSTOP® iX are covered by U.S. and international patents.

**LifeScience PLUS**  
REVOLUTIONARY WOUND CARE

Rev. 2019-10